



DENDIA

DENTAL INSTRUMENTS

...Dental Polishers

Perfect high gloss polishing of ZIRCONIUM OXIDE restorations


Perfekter Hochglanz für Zirkondioxid Restaurationen



DENDIA	ZIR9861M	ZIR9861F	ZIR9862M	ZIR9862F	ZIR9863M	ZIR9863F
Shank	RA	RA	RA	RA	RA	RA
Ø	040	040	060	060	140	140
L mm	10.0	10.0	7.5	7.5	1.6	1.6
Speed	7.000 - 10.000 rpm					
Packaging	2	2	2	2	2	2

2 - step diamond polishing system for removal, smoothing and gloss polishing

2 - Stufen Diamantpoliersystem zum Abtragen, Glätten und Hochglanzpolieren



Perfect high gloss polishing of COMPOSITE restorations

Perfekter Hochglanz für Komposit Restaurationen



DENDIA	9834	9832	9833	9837
Shank	RA	RA	RA	RA
Ø	040	060	100	140
L mm	10.0	9.5	0.7	1.6
Speed	7.000 - 10.000 rpm smoothing/glätten 3.000 - 8.000 rpm polishing/polieren			
Packaging	2	2	2	2


1 - step diamond polishing system for removal, smoothing and gloss polishing

1 - Stufen Diamantpoliersystem zum Abtragen, Glätten und Hochglanzpolieren



Perfect high gloss polishing of ALL-CERAMIC restorations

Perfekter Hochglanz für Vollkeramik Restaurationen



DENDIA	LS9871M	LS9871F	LS9872M	LS9872F	LS9873M	LS9873F
Shank	RA	RA	RA	RA	RA	RA
Ø	040	040	060	060	140	140
L mm	10.0	10.0	7.5	7.5	1.6	1.6
Speed	7.000 - 10.000 rpm					
Packaging	2	2	2	2	2	2

2 - step diamond polishing system for removal, smoothing and gloss polishing

2 - Stufen Diamantpoliersystem zum Abtragen, Glätten und Hochglanzpolieren



Silicon Carbide Polishers

Siliziumkarbid Polierer



DENDIA	9121M	9121F	9123M	9123F	9501G	9501M	9501F	9133G	9133M	9133F	9134G	9134M	9134F
Shank	RA	RA	RA	RA	RA	RA	RA	RA	RA	RA	RA	RA	RA
Ø	050	050	070	070	030	030	030	050	050	050	070	070	070
L mm	10.0	10.0	10.0	10.0	6.0	6.0	6.0	10.0	10.0	10.0	10.0	10.0	10.0
Speed	5.000 - 8.000 rpm				7.000 - 12.000 rpm			7.000 - 15.000 rpm					
Packaging	5				5			5					
Application	Composite				Amalgam/Gold			Ceramic					

The attractively priced alternative to diamond polishers
Die preislich attraktive Alternative zu Diamantpolierern

DENDIA

DENTAL INSTRUMENTS

Your local dealer:



CE 0483



DENDIA GmbH
Reichsstrasse 126
6800 Feldkirch
Austria
www.dendiadental.com